



BC32

Automated Blood Culture System


Even smaller laboratories can perform blood culture automatically because of me







- More Compact
- More Cost-Effective
- More Easily Controlled





3 Barcode scanning input

Colorimetric technology


Designed to detect microbial growth from blood specimens or sterile body fluids

Automated continuous monitoring for every 10 minutes for each bottle by independent detection sensor

Real-time display of incubation temperature and reaction chart

Flexible and quick historical data query and statistics functions





Four types of culture media available

Aerobic for Adult, Anaerobic for Adult, Aerobic for Pediatric and L-Form









Render Biotech Co., Ltd.

Add: 6K#, Building No.1, Baisha Emerging Industry Zone, Shahexi Road, Nanshan District, 518052, Shenzhen, China

Tel.: +86 755 8653 6648 E-mail: info@render-health.com

Web: www.renderbio.com